



# LoggerPro Video Analysis





- 1. Add a video.** Go to the Insert menu and select either type of video: a) “Movie...” to use an existing video. Browse to the desired file and select “Open” OR b) “Video Capture...” to record a new video using your camera. In the window, click “Start Video Capture” button to start recording and “Stop Video Capture” to stop.


---

- 2. Show tools.** Click **Enable analysis (A)**  on the lower right side of the window.

---

- 3. Sync up the time.** Use **Playback controls (B)**  allow you to play, stop, rewind, and increment next/previous frame. Navigate to the frame where the motion FIRST STARTS TO BEGIN. At this point, you can use the **Sync movie to graph (C)**  button and the “Graph time:” as 0 s. Leave the “Movie time” where it is.


---

- 4. Set the origin.** We want to set the release point as the origin ( $x=0$ ). When timed to release, use the **Set origin (D)**  tool to click a trackable part of the object. That sets that position as  $x=0$ .


---

- 5. Rotate the axes.** For an object moving only in one direction, we want that to be the x-direction. Depending on the camera angle, it may not be directly left and right. Click and rotate the yellow circle on the x-axis to **rotate the axes** so they align with the path of motion.


---


- 6. Set the scale.** Using **Set scale (E)** , click and drag along a known length of an object. Enter in the box the actual size of the object used.

---

- 7. Add data points.** Find a part of the cart that is easily visible and trackable. Using the **Add points (F)**  tool, click on that point to **add a data point** and advance the video by one frame. **Repeat clicking until data collection is complete.** (When object is touched to stop or reaches the end).

---

- 8. Add best fit lines.** In the window’s shortcut bar, use the **Curve fit tool**  to add a best fit line (linear, quadratic, inverse, etc). To analyze motion in the x-direction, when prompted select “VideoAnalysis|X”. Select the desired model in the window under “general equation”. Click “test fit” to get a quick visual to make sure the model fits the data, then select “OK.”

If you know the fit is linear, you can use the **linear fit tool**  to save some steps.